


UVa/Padova
T1DMS
Type 1 Diabetes Metabolic Simulator
Installation Guide v3.2


Are you upgrading your from T1DMS v2? STOP. Launch the License Manager. De-Activate your existing license first! Then Activate your License. Your License # will now work for both v2 and v3.2.

Please follow the instructions below to install the UVa/Padova T1DM Simulator v3.2:


1. Make sure that you exit all other applications before installation
2. Check with your IT Personnel to make sure you can install software on your hard drive
3. Double-click on the installer **UVa_Padova_T1DM_Simulator_Setup_v3.2.exe**:




4. Click **Next** if you wish to install T1DMS:




5. Click **Next** to confirm and proceed with installation:




6. Please read the terms and conditions, then select **I agree with the above terms and conditions** then click **Next** to proceed:




7. Enter the name of the installation folder on your hard drive then click **Next**:



8. If the folder does not exist, you may create the folder by clicking **Yes**:




9. Click **Start** to initiate the installation process:




Are you upgrading your license from v2? STOP. When you Launch the License Manager. De-Activate your existing license first! Then Activate your License. Your License # will now work for both v2 and v3.2.

10. Click **Exit** to finish the installation and launch the License Manager to activate your software:




Please follow the instructions below to activate the T1DMS software:

11. Click **Yes** to run the License Manager to activate your license and register your software:




12. Enter your 18-digit License Number *exactly as provided* and click on **Activate**:

(Note: if you do not have your License Number, please close this window and activate your license at a later date by going to your Program Files Folder for the T1DMS Simulator software and click on your License Manager to return to the activation window)



13. If you have not previously registered the T1DMS software, click **No** to create an account. If you have previously registered, click **Yes** to sign in:

(Note: you must create a username and password to allow for future access to your account)





14. If you clicked **No** above to create an account, you will be directed to the Create Account window to provide account information then click on **OK**:



(Note: the username and password created here will be needed for future re-activation if the license is de-activated)

Please fill out the following information. Fields in bold are required.

Display Language	en - English UK
First Name	
Last Name	
Address 1	
Address 2	
Address 3	
City/Town	
State/Province	
Zip/Postal Code	
Phone 1	
Phone 2	
Fax 1	
Fax 2	
Country	Abkhazia
Company Name	
Position in Company	
Username	
E-mail address	
	<input type="checkbox"/> Yes add me to the email database
	<input type="checkbox"/> Yes I want to receive offers from selected partners
E-mail Format	both
PASSWORD	
Confirm Password	
	OK

*(Special Note: Your computer must have direct access to the internet to activate & de-activate T1DMS Simulator.
Import Certificate/Export Certificate activation is not supported)*


15. Once the software has been activated, click on **OK**:




Please follow the instructions below to de-activate your license:

You can de-activate your T1DMS license at any time. This action returns your license to the server for storage and re-use, allowing you to quickly and easily move the license to another machine. Once de-activated, the T1DMS software will no longer run on the machine until you (1) re-activate it on the same machine, or (2) install the T1DMS software on a new machine and activate it on that machine. This deactivation/activation process allows you to de-activate and activate between two computers easily:


1. Go to your Program Files Folder for the T1DMS software and double-click on **License Manager.exe**:




2. Once the License Manager window opens, click on **De-activate** in the Menu Bar:



3. Enter your 18-digit License code and click on **De-Activate**:




4. Your license can now be activated on another computer or re-activated on the same computer. Click on **OK** to proceed:



Please follow the instructions below to re-activate your license:


1. If you try to run the UVa/Padova T1DM Simulator software without re-activating the license, you will receive the following message:




and a Matlab Error:




2. Go to your Program Files Folder for the T1DMS software and double-click on **License Manager.exe**:




3. Enter your 18-digit License code and click **Activate**:





4. Sign into your account by clicking on **Yes**:



5. Enter your Username and Password then click on **OK**:



6. You will receive the following two messages – click **OK** on both to complete re-activation:





The Epsilon Group
2000 Holiday Drive
Charlottesville, VA 22901 USA

Telephone: (434) 923-6091

info@tegvirginia.com
www.tegvirginia.com

Follow TEG on:



<http://www.linkedin.com/company/2653952>



<http://www.twitter.com/tegvirginia>