

Guide to Installing T1DMS Software

The Epsilon Group, an ARDx Informatics
Company
2000 Holiday Drive
Charlottesville, Va 22901

<https://tegvirginia.com/contact/>



© 2020 ARDx Informatics. All rights reserved.

Version Identification

Version	Date	Description	Author
1	07/30/2020	Installer for T1DMS version 3.2.1.	K. Bowden


UVa/Padova T1DMS

Type 1 Diabetes Metabolic Simulator


Installation Guide v3.2

Please follow the instructions below to install the UVa/Padova T1DM Simulator v3.2:


1. Make sure that you exit all other applications before installation
2. Check with your IT Personnel to make sure you can install software on your hard drive.
3. Double-click on the installer **UVa_Padova_T1DM_Simulator_Setup_v3.2.exe**:




4. Click **Next** if you wish to install T1DMS:




5. Click **Next** to confirm and proceed with installation:




6. Please read the terms and conditions, then select **I agree with the above terms and conditions** then click **Next** to proceed:




7. Enter the name of the installation folder on your hard drive then click **Next**:



8. You are now ready to install UVA/PADOVA T1DM Simulator v3.2.1. Click Start:




9. Click **Exit** to finish the installation and launch the License Manager to activate your software:



Please follow the instructions below to activate the T1DMS software:

10. This is the License Manager Screen to activate your software:



TEG License Manager

Initializing license system...
License system initialized successfully.
License status is shown in the "About" tab.

Activation Deactivation About

☒ Online ☐ Offline Enter activation code and required customer information for online activation (Internet connection is required)

Activation Code

Customer Data (* is required)


Name* E-mail*

Phone Company

Activate

Note: If your computer does not have access to the internet, you will need to use the offline activation process. If this applies to you, please contact The Epsilon Group (info@tegvirginia.com) for assistance.

11. Enter the 16-digit License Number, Your Name and E-mail address. Click on **Activate**:



TEG License Manager

Initializing license system...
License system initialized successfully.
License status is shown in the "About" tab.
Ready for online activation.
Ready for online activation.

Activation Deactivation About

☒ Online ☐ Offline Enter activation code and required customer information for online activation (Internet connection is required)

Activation Code


Customer Data (* is required)

Name* E-mail*


Phone Company

Activate

12. Your license is now activated. This screen will notify you (at the top) of the activation.




13. Go to your Program Files Folder for the T1DMS software and double-click on Run Simulator:




Please follow the instructions below to de-activate your license:

You can de-activate your T1DMS license at any time. This action returns your license to the server for storage and re-use, allowing you to quickly and easily move the license to another machine. Once de-activated, the T1DMS software will no longer be available until you (1) re-activate it on the same machine, or (2) install the T1DMS software on a new machine and activate it on that machine. This deactivation/activation process allows you to de-activate and activate between two computers easily:


1. Go to your Program Files Folder for the T1DMS software and double-click on **License Manager**:



2. Once the License Manager window opens, click on **De-activation** in the Menu Bar:




3. Your License has now been **de-activated** (“revoked” as shown in messages at top of screen):




Please follow the instructions below to re-activate your license:

1. Go to your Program Files Folder for the T1DMS software and double-click on **License Manager**:



2. This is the License Manager Screen to re-activate your software:



TEG License Manager

Initializing license system...
License system initialized successfully.
License status is shown in the "About" tab.

Activation **Deactivation** **About**

☒ Online ☐ Offline Enter activation code and required customer information for online activation (Internet connection is required)


Activation Code

Customer Data (* is required)

Name* E-mail*

Phone Company

3. Enter the 16-digit License Number, Your Name and E-mail address. Click on Activation:



TEG License Manager

Initializing license system...
License system initialized successfully.
License status is shown in the "About" tab.
Ready for online activation.
Ready for online activation.

Activation **Deactivation** **About**

☒ Online ☐ Offline Enter activation code and required customer information for online activation (Internet connection is required)

Activation Code: 4068-6391-2954-5596


Customer Data (* is required)

Name* Test Person C E-mail* kbowden@tegvirginia.com

Phone Company

Activate

4. Your license is now activated. This screen will notify you (at top of screen) of the activation.



TEG License Manager

Initializing license system...
License system initialized successfully.
License status is shown in the "About" tab.
Ready for online activation.
Ready for online activation.
Activating, please wait...
Ready to activate with activation key: 4068-6391-2954-5596
Validating license file...
Write license file succeeded.
Internet activation completed.
The activation code was activated successfully.
Validating the license file...
Activation succeeded.
Congratulations! T1DMS 3.2.1 is licensed.

Activation **Deactivation** **About**

☒ Online ☐ Offline Enter activation code and required customer information for online activation (Internet connection is required)

Activation Code: 4068-6391-2954-5596

Customer Data (* is required)

Name* Test Person C E-mail* kbowden@tegvirginia.com

Phone Company

Activate

**The Epsilon Group,
an Abbott Rapid Diagnostics Informatics Company**
2000 Holiday Drive
Charlottesville, VA 22901 USA

Telephone: (434) 923-6091

E-mail: info@tegvirginia.com
www.tegvirginia.com